

Agrobacterium infiltration of *Nicotiana benthamiana* for transient protein expression – Catherine Espinoza Patharkar

Transient expression in the wild tobacco plant (*Nicotiana benthamiana*) is used to determine the subcellular location of a protein of interest when tagged with a reporter such as green fluorescent protein (GFP), or to mass produce proteins without making transgenic plants. A domesticated version of the crown gall forming bacteria, *Agrobacterium tumefaciens*, is used to introduce the target gene expression cassette into *N. benthamiana* leaves.

Materials

- *Agrobacterium* strain hosting a plant expression construct (usually driven by Cauliflower mosaic virus 35S promoter)
- Healthy non-flowering tobacco (*Nicotiana benthamiana*) plants (3-4 weeks old)
- MES/KOH (pH 5.6)
- MgCl₂
- Acetosyringone
- LB media with appropriate antibiotics

The time frame for the protocol described here is as follows:

- Day 1: *Agrobacterium* transformation with binary vectors.
- Day 3: Start *Agrobacterium* culture (including *Agrobacterium* containing p19 strain). *p19* is a suppressor of gene silencing from the tomato bushy stunt virus and must be coinfiltrated for proper gene expression. Water tobacco plants.
- Day 4: Infiltration of transformed *Agrobacteria* in *N. benthamiana* leaves
- Day 6-10: Imaging of results.

Agrobacterium*-Mediated Infiltration of *N. benthamiana

This method was modified from [Voignet et al. \(2003\)](#).

1. Prepare the Activation buffer

Reagent stock concentration	10 ml	Final concentration
1 M MES/KOH (pH 5.6)	100 μ l	10 mM
10 mM MgCl ₂	10 ml	10 mM
150 mM Acetosyringone (3',5'-Dimethoxy-4'-hydroxyacetophenone) in DMSO, stored at -20°C	10 μ l	150 μ M

2. Measure the OD₆₀₀ of the overnight cultures.
3. *N. benthamiana* will be infiltrated with a solution containing OD₆₀₀=0.3 of each experimental Agro containing construct and OD₆₀₀=0.1 of p19. Calculate the volume of cultures ($V_{\text{construct}}$; V_{p19}) needed according to the formulas:

$V_{\text{construct}} = V_{\text{final}} \times 0.3/\text{OD}_{600}$; $V_{\text{p19}} = V_{\text{final}} \times 0.1/\text{OD}_{600}$. One mL of infiltrate is often enough to complete a small experiment. Plan your final volume accordingly.

4. Mix calculated volumes of construct(s) and p19 cultures into a 1.5 ml Eppendorf tube.
5. Immediately centrifuge the mixtures at maximum speed for 30 sec at room temperature.
6. Discard the supernatant and resuspend the mixtures in V_{final} of activation buffer.
7. Incubate the mixtures for at least 30 min at room temperature in dark.
8. Infiltrate the mixtures using a 1-mL syringe without a needle into the abaxial (bottom) side of *N. benthamiana* leaves. Trace the infiltrated area with a marker. Incubate the plants in the greenhouse for 2-6 d, depending on the level of protein expression. Two days of expression is often maximal.

A video protocol can be viewed at: <https://rahulpatharkar.000webhostapp.com/?p=313>